

Prescription for Home Delivery



Please complete this form and return via fax (855-699-3366) or email (info@pawprintoxygen.com).
If this form is received M-F from 9am-4pm EST, we will contact the pet owner within 1 hour.

PATIENT INFORMATION

Owner Name: _____ Patient Name: _____ Age: _____
Owner Phone: _____ Species/Breed: _____
Owner Email: _____ Weight: _____ Sex: _____
Owner Address: _____ Relevant Conditions: _____

PRESCRIPTION INFORMATION

Portable Oxygen Kit

Portable Oxygen (99%), Flow by Mask

Choose a Kit:

- Dog, 0-11 lb (0-5 kg) Cat, 0-11 lb (0-5 kg)
 Dog, 11-22 lb (6-10 kg) Cat, 11+ lb (6+ kg)
 Dog, 22-33 lb (11-15 kg) Custom Kit
 Dog, 33+ lb (16+ kg) *Write flow rate below, see ➤ section on back*

Instructions for use: _____

Oxygen Concentrator

Stationary Oxygen (>95%) at 1-10 LPM, Flow by Mask or Cage

via Oxygen Cage - Indicate desired oxygen level:

- 25-35% 35-45% 45-55%

See Page 2 for details on the Oxygen Cage and how the venturi system is used to control the oxygen level in the Cage.

via Pet Oxygen Mask - Indicate flow rate: _____

Instructions for use: _____

HOSPITAL INFORMATION

Hospital Name: _____ Contact Number: _____
License #: _____ State: _____ Contact Email: _____
Authorizing Doctor: _____ (print) _____ (sign) _____ (date)



Questions about this form?
Scan the QR code or call 855-699-4366

Phone: (855) 699-4366
Fax: (855) 699-3366
info@pawprintoxygen.com

PRESCRIPTION GUIDELINES

Veterinarians: You are responsible for determining the oxygen flow rate and/or oxygen level based on the delivery method, animal weight, and desired inspired concentration. While Pawprint Oxygen provides guidelines, it is the attending veterinarian's responsibility to determine administration and dosage.

Portable Oxygen Kit

A portable supply of oxygen for use at home and during transport. Each Kit includes Pawprint Oxygen Canisters, Regulator, PureVent Mask and Tubing.

Standard Kits:

Kit	Mask	Flow Rate	Duration
Cat 0-11 lb	Small	0.5 LPM	90 min
Cat 11+ lb	Small	1.0 LPM	90 min
Dog 0-11 lb	Medium	0.5 LPM	90 min
Dog 11-22 lb	Medium	1.0 LPM	90 min
Dog 22-33 lb	Large	2.0 LPM	45 min
Dog 33+ lb	Large	3.0 LPM	30 min

The **Flow Rate** for standard kits is based on a dosage of 100 ml/kg/min. The Regulator is factory-set to the prescribed flow rate and cannot be changed by the client.

The **Duration** of the kit can be increased by purchasing additional Pawprint Oxygen Canisters, which contain 15 liters of oxygen (99%) each.

The **Mask** refers to the PureVent Pet Oxygen Mask, which is included in all kits.

➤ **Custom Kits:** Write in the flow rate (0.5, 1.0, 2.0, or 3.0 LPM). The maximum flow rate is 3.0 LPM. If you wish to prescribe a higher flow rate, please consider an Oxygen Concentrator. For custom kits, the correct Mask size will be provided by our experts based on the pet.

Scan the QR code below for more information, including dosage guidelines, how to prescribe custom kits for exotics and other animals, and more FAQs.

You can order Pawprint Oxygen products for your practice from MWI Animal Health, Vetcove, or by calling us directly at 855-699-4366.



Questions about this form?
Scan the QR code or call 855-699-4366



Oxygen Concentrator

Oxygen Concentrators supply a continuous flow of medical-grade oxygen (>95%) and may be used with an Oxygen Cage and/or a Pet Oxygen Mask.

Administration via Oxygen Cage

An Oxygen Cage may be prescribed for use with the Oxygen Concentrator. The oxygen level in the Cage can be set between 25-55% using the included venturi pieces. **To approve the use of an Oxygen Cage, veterinarians must indicate the desired oxygen level.** The actual oxygen level may vary by +/-5 % due to venting and air exchange.

Oxygen Cages are available in the following sizes:

- Medium (24" x 18" x 18"): For pets under 25 lbs
- Large (44" x 26" x 26"): For pets under 65 lbs

Administration via Pet Oxygen Mask

A PureVent Pet Oxygen Mask may be purchased with the Oxygen Concentrator. Vets should indicate a flow rate (1-10 LPM) for administration by Mask.

Ordering Process

Pawprint Oxygen is open Monday-Friday 9am-5pm EST. If you are submitting this form within normal business hours, pet owners can expect a call from our team within 1 hour to complete their order. If you are submitting this form outside of normal business hours, pet owners can expect a call the next business day.

Pawprint Oxygen handles all billing, shipping, and equipment education directly with the client for home delivery orders.

If this request is urgent or requires special instructions, please call us at 855-699-4366.

Phone: (855) 699-4366

Fax: (855) 699-3366

info@pawprintoxygen.com