

Prescription Approval Form

Please complete Page 1 of this form and return it via fax (855-699-3366) or email (info@pawprintoxygen.com). If this form is received Monday-Friday between 9am-5pm EST, we will contact the pet owner within 2 hours.

Pawprint Oxygen Canisters

Oxygen USP (99%), Flow by Mask

PLEASE INDICATE APPROVAL:

APPROVED NOT APPROVED

PLEASE INDICATE FLOW RATE:

0.5 LPM 1.0 LPM 2.0 LPM

Includes portable Canisters (10 liters of oxygen per Canister), a Regulator, and a Pet Oxygen Mask. Best for episodic use or for transport. See Page 2 for details on Pawprint Oxygen Canisters and Mask sizing.

Extended Oxygen Therapy

Oxygen (>95%), Flow by Mask or Cage

PLEASE INDICATE APPROVAL:

APPROVED NOT APPROVED

Includes an Oxygen Concentrator (>95% purity @ 1-10 LPM), an Oxygen Cage, and a Pet Oxygen Mask. Flow rate depends on the delivery method. Pet owners will be informed using guidelines on Page 2. Veterinarians may write further instructions below:

Patient Information

Owner Name: _____

Patient Name: _____ Age: _____

Owner Phone: _____

Species/Breed: _____

Owner Email: _____

Weight: _____ Gender: _____

Owner Address: _____

Relevant Conditions: _____

Hospital Information

Hospital Name: _____

Contact Number: _____

License #: _____ State: _____

Contact Email: _____

Authorizing Doctor: _____

(print)

(sign)

(date)



Prescription Guidelines

Veterinarians: You are responsible for determining the oxygen flow rate based on the delivery method, animal weight, and desired inspired concentration. While Pawprint Oxygen provides guidelines (see chart below), it is the attending veterinarian's responsibility to determine dosage and oxygen flow rate.

Pawprint Oxygen Canisters

Pawprint Oxygen Canisters are recommended for pets that need relief from episodic respiratory distress. This is a portable option that can be used for on-the-spot treatment at home and on the go.

Pawprint Oxygen Canisters are most commonly used in tandem with the Pet Oxygen Mask (included). The Mask is available in Small, Medium, and Large sizes. We will provide the correct size based on the patient.

Three flow rates are available: 0.5, 1.0, and 2.0 LPM. The Regulator will be factory set to the flow rate prescribed by the veterinarian and **is not** adjustable once it is set. If no flow rate is indicated, the reference chart below will be used to determine the flow rate:

Dosage: 50 mL/kg/min (nasal cannula)

Animal Weight	Flow Rate	Time Till Empty
0-10 kg (0-22 lb)	0.5 LPM	20 min
11-20 kg (23-44 lb)	1.0 LPM	10 min
21-40 kg (45-88 lb)	2.0 LPM	5 min

Dosage: 100 mL/kg/min (Pet Oxygen Mask)

Animal Weight	Flow Rate	Time Till Empty
0-5 kg (0-11 lb)	0.5 LPM	20 min
6-10 kg (12-22 lb)	1.0 LPM	10 min
11-15 kg (23-33 lb)	2.0 LPM	5 min

Dosage: 200 mL/kg/min (flow by)

Animal Weight	Flow Rate	Time Till Empty
0-5 kg (0-11 lb)	1.0 LPM	10 min
6-10 kg (12-22 lb)	2.0 LPM	5 min

Extended Oxygen Therapy

Oxygen Concentrators may be used for pets that require more than one hour of therapy per day or pets requiring higher flow rates than can be achieved with Oxygen Canisters. Concentrators may be used in tandem with a Mask or a Cage.

The Concentrator is available in two models:

- 5 Liter: Adjustable flow rate from 1-5 LPM
- 10 Liter: Adjustable flow rate from 2-10 LPM

Oxygen Cages are available in the following sizes:

- Small (18" x 13" x 13"): For pets under 5 lbs
- Medium (24" x 18" x 18"): For pets under 25 lbs
- Large (44" x 26" x 26"): For pets under 75 lbs

To achieve sufficient oxygen saturation levels:

- Use a 5 Liter Concentrator to supply oxygen to a Small or Medium Cage.
- Use a 10 Liter Concentrator to supply oxygen to a Large Cage.

Note: Pawprint Oxygen Canisters **cannot** be used to supply Oxygen Cages due to flow rate requirements.

Ordering Process

Pawprint Oxygen is open Monday-Friday 9am-5pm EST. If you are submitting this form within normal business hours, pet owners can expect a call from our team within 2 hours to complete their order. If you are submitting this form outside of normal business hours, pet owners can expect a call the next business day.

If this is urgent or requires special instructions, please call us at 855-699-4366. If you would like to dispense Pawprint Oxygen directly from your practice, please call or visit pawprintoxygen.com

